

ABX Minocal

LOT

CX 465

CAL



(Exp.) 2022-01-05

Rev 1

(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio								TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA		PENTRA
		ES60	CRP 200			60 8P 60 18P					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.9	9.0			9.0					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.43	4.44			4.38					± 0.06
HB HGB	g/dl	13.6	13.4			13.6					± 0.2
	g/l	136	134			136					± 2
	mmol/l	8.45	8.34			8.45					± 0.12
HT HCT	%	38.5	38.6			39.4					± 1.0
	l/l	0.385	0.386			0.394					± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	244	250			256					± 10
VPM MPV	µm ³ ; fl	9.3	8.9			9.7					± 0.5

ABX Minocal

LOT

CX 465

CAL



(Exp.) 2022-01-05
(YYYY-MM-DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	LyseVet								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC	10 ⁹ /mm ³ ;10 ⁹ /l	8.5	8.6							± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.40	4.36							± 0.06
HB HGB	g/dl	13.2	13.2							± 0.2
	g/l	132	132							± 2
	mmol/l	8.20	8.20							± 0.12
HT HCT	%	37.8	37.5							± 1.0
	l/l	0.378	0.375							± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	231	245							± 10
VPM MPV	µm ³ ; fl	N/A	9.5							± 0.5

Ref: TEMP-0387 Rev.47 FRONT / RECTO 1300060881