

ABX Minocal

LOT

CX 464

CAL



(Exp.) 2021-12-05

Rev 1

(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio								TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA		PENTRA
		ES60	CRP 200			60 8P 60 18P					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.9	9.2			9.0					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.49	4.55			4.47					± 0.06
HB HGB	g/dl	13.4	13.3			13.5					± 0.2
	g/l	134	133			135					± 2
	mmol/l	8.32	8.28			8.38					± 0.12
HT HCT	%	37.7	38.7			39.3					± 1.0
	l/l	0.377	0.387			0.393					± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	237	248			244					± 10
VPM MPV	µm ³ ; fl	9.0	8.6			9.3					± 0.5

ABX Minocal

LOT

CX 464

CAL



(Exp.) 2021-12-05
(YYYY-MM-DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	LyseVet								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC	10 ⁹ /mm ³ ;10 ⁹ /l	8.7	8.6							± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.48	4.47							± 0.06
HB HGB	g/dl	13.1	13.1							± 0.2
	g/l	131	131							± 2
	mmol/l	8.14	8.14							± 0.12
HT HCT	%	37.6	37.5							± 1.0
	l/l	0.376	0.375							± 0.010
PLAQ. PLTS	10 ⁹ /mm ³ ;10 ⁹ /l	235	243							± 10
VPM MPV	µm ³ ; fl	N/A	9.4							± 0.5

Ref: TEMP-0387 Rev.47 FRONT / RECTO 1300060881