

# ABX Minocal

**LOT**

**CX 464**

**CAL**



**(Exp.) 2021-12-05**

Rev 1

(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio								TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA		PENTRA
		ES60	CRP 200			60 8P 60 18P					
GB WBC	10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	<b>8.9</b>	<b>9.2</b>			<b>9.0</b>					± 0.2
GR RBC	10 <sup>6</sup> /mm <sup>3</sup> ;10 <sup>12</sup> /l	<b>4.49</b>	<b>4.55</b>			<b>4.47</b>					± 0.06
HB HGB	g/dl	<b>13.4</b>	<b>13.3</b>			<b>13.5</b>					± 0.2
	g/l	<b>134</b>	<b>133</b>			<b>135</b>					± 2
	mmol/l	<b>8.32</b>	<b>8.28</b>			<b>8.38</b>					± 0.12
HT HCT	%	<b>37.7</b>	<b>38.7</b>			<b>39.3</b>					± 1.0
	l/l	<b>0.377</b>	<b>0.387</b>			<b>0.393</b>					± 0.010
PLAQ. PLTS	10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	<b>237</b>	<b>248</b>			<b>244</b>					± 10
VPM MPV	µm <sup>3</sup> ; fl	<b>9.0</b>	<b>8.6</b>			<b>9.3</b>					± 0.5

# ABX Minocal

**LOT**

**CX 464**

**CAL**



**(Exp.) 2021-12-05**  
(YYYY-MM-DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	LyseVet								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC	10 <sup>9</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	<b>8.7</b>	<b>8.6</b>							± 0.2
GR RBC	10 <sup>6</sup> /mm <sup>3</sup> ;10 <sup>12</sup> /l	<b>4.48</b>	<b>4.47</b>							± 0.06
HB HGB	g/dl	<b>13.1</b>	<b>13.1</b>							± 0.2
	g/l	<b>131</b>	<b>131</b>							± 2
	mmol/l	<b>8.14</b>	<b>8.14</b>							± 0.12
HT HCT	%	<b>37.6</b>	<b>37.5</b>							± 1.0
	l/l	<b>0.376</b>	<b>0.375</b>							± 0.010
PLAQ. PLTS	10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	<b>235</b>	<b>243</b>							± 10
VPM MPV	µm <sup>3</sup> ; fl	N/A	<b>9.4</b>							± 0.5

Ref: TEMP-0387 Rev.47 FRONT / RECTO 1300060881