

ABX Minocal

LOT

CX 419

CAL



(Exp.) 2018-03-05

Rev 1

(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio								TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA		PENTRA
		ES60	CRP 200			60 8P 60 18P					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.1	9.4			9.4					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.52	4.54			4.51					± 0.06
HB HGB	g/dl	13.9	13.6			13.8					± 0.2
	g/l	139	136			138					± 2
	mmol/l	8.63	8.46			8.57					± 0.12
HT HCT	%	39.2	39.5			39.2					± 1.0
	l/l	0.392	0.395			0.392					± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	255	250			258					± 10
VPM MPV	µm ³ ; fl	9.2	8.8			9.0					± 0.5

ABX Minocal

LOT

CX 419

CAL



(Exp.) 2018-03-05
(YYYY-MM-DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	LyseVet								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC	10 ⁹ /mm ³ ;10 ⁹ /l	8.9	9.2							± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.53	4.53							± 0.06
HB HGB	g/dl	13.6	13.5							± 0.2
	g/l	136	135							± 2
	mmol/l	8.45	8.38							± 0.12
HT HCT	%	39.2	39.2							± 1.0
	l/l	0.392	0.392							± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	235	245							± 10
VPM MPV	µm ³ ; fl	N/A	9.0							± 0.5

Ref: TEMP-0387 Rev.44 FRONT / RECTO 1300032436-A