

# ABX Minocal



**LOT** CX 071  
Rev 1

**CAL**

(Exp.) 2011-08-05  
(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX VETPACK										TOLERANCES TOLERANCE
		MICROS	MICROS									
		Scil Vet ABC	Scil Vet ABC Plus ESV 60									
GB WBC	10 <sup>3</sup> /mm <sup>3</sup> ; 10 <sup>9</sup> /l	<b>9.2</b>	<b>9.3</b>									± 0.2
GR RBC	10 <sup>6</sup> /mm <sup>3</sup> ; 10 <sup>12</sup> /l	<b>4.63</b>	<b>4.71</b>									± 0.06
HB HGB	g/dl	<b>13.5</b>	<b>13.3</b>									± 0.2
	g/l	<b>135</b>	<b>133</b>									± 2
	mmol/l	<b>8.38</b>	<b>8.28</b>									± 0.12
HT HCT	%	<b>39.4</b>	<b>39.6</b>									± 1.0
	l/l	<b>0.394</b>	<b>0.396</b>									± 0.010
PLAQ. PLTS	10 <sup>9</sup> /mm <sup>3</sup> ; 10 <sup>9</sup> /l	<b>243</b>	<b>239</b>									± 10
VPM MPV	µm <sup>3</sup> ; fl	N/A	<b>8.6</b>									± 0.5

Ref: FORM-0387 Rev.19 FRONT / RECTO 9930044-B