

# ABX Minocal



**LOT** CX 123  
Rev 1

**CAL**

 (Exp.) 2014-01-05  
(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX VETPACK										TOLERANCES TOLERANCE
		MICROS	MICROS									
		Scil Vet ABC	Scil Vet ABC Plus ESV 60									
GB WBC	10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	8.9	9.0									± 0.2
GR RBC	10 <sup>6</sup> /mm <sup>3</sup> ;10 <sup>12</sup> /l	4.47	4.50									± 0.06
HB HGB	g/dl	13.4	13.4									± 0.2
	g/l	134	134									± 2
	mmol/l	8.32	8.32									± 0.12
HT HCT	%	38.0	38.5									± 1.0
	l/l	0.380	0.385									± 0.010
PLAQ. PLTS	10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	238	243									± 10
VPM MPV	µm <sup>3</sup> ; fl	N/A	8.1									± 0.5

Ref: TEMP-0387 Rev.33 9930078-A