

ABX Minocal



LOT CX 103
Rev 1

CAL

 (Exp.) 2013-11-05
(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX VETPACK										TOLERANCES TOLERANCE
		MICROS	MICROS									
		Scil Vet ABC	Scil Vet ABC Plus ESV 60									
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.1	9.2									± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.62	4.59									± 0.06
HB HGB	g/dl	13.5	13.3									± 0.2
	g/l	135	133									± 2
	mmol/l	8.38	8.26									± 0.12
HT HCT	%	38.8	38.6									± 1.0
	l/l	0.388	0.386									± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	242	238									± 10
VPM MPV	µm ³ ; fl	N/A	8.3									± 0.5

Ref: TEMP-0387 Rev.33 9930078-A