

# ABX Minocal



**LOT** CX 093  
Rev 1

**CAL**

(Exp.) 2013-10-05  
(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX VETPACK										TOLERANCES TOLERANCE
		MICROS	MICROS									
		Scil Vet ABC	Scil Vet ABC Plus ESV 60									
GB WBC	10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	8.9	8.9									± 0.2
GR RBC	10 <sup>6</sup> /mm <sup>3</sup> ;10 <sup>12</sup> /l	4.51	4.47									± 0.06
HB HGB	g/dl	13.5	13.3									± 0.2
	g/l	135	133									± 2
	mmol/l	8.38	8.26									± 0.12
HT HCT	%	38.8	38.4									± 1.0
	l/l	0.388	0.384									± 0.010
PLAQ. PLTS	10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	247	246									± 10
VPM MPV	µm <sup>3</sup> ; fl	N/A	8.5									± 0.5

Ref: TEMP-0387 Rev.33 9930078-A