

ABX Minocal



LOT CX 083
Rev 1

CAL

(Exp.) 2013-09-05
(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX VETPACK										TOLERANCES TOLERANCE
		MICROS		MICROS								
		Scil Vet ABC	Scil Vet ABC Plus ESV 60									
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.6	8.7									± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.52	4.52									± 0.06
HB HGB	g/dl	13.4	13.3									± 0.2
	g/l	134	133									± 2
	mmol/l	8.32	8.26									± 0.12
HT HCT	%	38.6	38.7									± 1.0
	l/l	0.386	0.387									± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	230	235									± 10
VPM MPV	µm ³ ; fl	N/A	8.7									± 0.5

Ref: TEMP-0387 Rev.33 9930078-A