

ABX Minocal



LOT

CX 073

CAL



(Exp.) 2013-08-05

(YYYY - MM - DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	ABX VETPACK										TOLERANCES TOLERANCE
		MICROS	MICROS									
		Scil Vet ABC	Scil Vet ABC Plus ESV 60									
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.0	9.1									± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.45	4.46									± 0.06
HB HGB	g/dl	13.4	13.2									± 0.2
	g/l	134	132									± 2
	mmol/l	8.32	8.20									± 0.12
HT HCT	%	38.3	38.4									± 1.0
	l/l	0.383	0.384									± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	231	235									± 10
VPM MPV	µm ³ ; fl	N/A	8.5									± 0.5

Ref: TEMP-0387 Rev.33 9930078-A