

ABX Minocal



LOT CX 063
Rev 1

CAL

(Exp.) 2013-07-05
(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX VETPACK										TOLERANCES TOLERANCE
		MICROS	MICROS									
		Scil Vet ABC	Scil Vet ABC Plus ESV 60									
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.8	8.9									± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.50	4.48									± 0.06
HB HGB	g/dl	13.4	13.4									± 0.2
	g/l	134	134									± 2
	mmol/l	8.32	8.32									± 0.12
HT HCT	%	38.7	39.0									± 1.0
	l/l	0.387	0.390									± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	252	255									± 10
VPM MPV	µm ³ ; fl	N/A	8.4									± 0.5

Ref: TEMP-0387 Rev.33 9930078-A