

HORIBA ABX SAS
HORIBA ABX SAS
Parc Euromédecine
Rue du Caducée - BP 7290
34184 MONTPELLIER Cedex 4
FRANCE
Tél.: 33 (0) 4 67 14 15 16
Fax.: 33 (0) 4 67 14 15 17

LOT CX 053
Rev 1

CAL

 (Exp.) 2013-06-05
(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX VETPACK										TOLERANCES TOLERANCE
		MICROS	MICROS									
		Scil Vet ABC	Scil Vet ABC Plus ESV 60									
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.0	9.0									± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.54	4.50									± 0.06
HB HGB	g/dl	13.4	13.3									± 0.2
	g/l	134	133									± 2
	mmol/l	8.32	8.26									± 0.12
HT HCT	%	38.4	38.3									± 1.0
	l/l	0.384	0.383									± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	240	240									± 10
VPM MPV	µm ³ ; fl	N/A	7.8									± 0.5

Ref: TEMP-0387 Rev.32 9930078-A