

HORIBA ABX SAS
HORIBA ABX SAS
Parc Euromédecine
Rue du Caducée - BP 7290
34184 MONTPELLIER Cedex 4
FRANCE
Tél.: 33 (0) 4 67 14 15 16
Fax.: 33 (0) 4 67 14 15 17

LOT CX 043
Rev 1

CAL

(Exp.) 2013-05-05
(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX VETPACK										TOLERANCES TOLERANCE
		MICROS	MICROS									
		Scil Vet ABC	Scil Vet ABC Plus ESV 60									
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.3	9.3									± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.53	4.55									± 0.06
HB HGB	g/dl	13.5	13.4									± 0.2
	g/l	135	134									± 2
	mmol/l	8.38	8.32									± 0.12
HT HCT	%	38.9	38.9									± 1.0
	l/l	0.389	0.389									± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	243	247									± 10
VPM MPV	µm ³ ; fl	N/A	8.1									± 0.5

Ref: TEMP-0387 Rev.31 9930078-A