

**HORIBA ABX SAS**  
HORIBA ABX SAS  
Parc Euromédecine  
Rue du Caducée - BP 7290  
34184 MONTPELLIER Cedex 4  
FRANCE  
Tél.: 33 (0) 4 67 14 15 16  
Fax.: 33 (0) 4 67 14 15 17

**LOT** CX 033  
Rev 1

**CAL**

(Exp.) 2013-04-05  
(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX VETPACK										TOLERANCES TOLERANCE
		MICROS	MICROS									
		Scil Vet ABC	Scil Vet ABC Plus ESV 60									
GB WBC	10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	9.0	9.0									± 0.2
GR RBC	10 <sup>6</sup> /mm <sup>3</sup> ;10 <sup>12</sup> /l	4.63	4.60									± 0.06
HB HGB	g/dl	13.6	13.4									± 0.2
	g/l	136	134									± 2
	mmol/l	8.45	8.32									± 0.12
HT HCT	%	38.9	38.9									± 1.0
	l/l	0.389	0.389									± 0.010
PLAQ. PLTS	10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	234	235									± 10
VPM MPV	µm <sup>3</sup> ; fl	N/A	8.1									± 0.5

Ref: TEMP-0387 Rev.31 9930078-A