

HORIBA ABX SAS
HORIBA ABX SAS
Parc Euromédecine
Rue du Caducée - BP 7290
34184 MONTPELLIER Cedex 4
FRANCE
Tél.: 33 (0) 4 67 14 15 16
Fax.: 33 (0) 4 67 14 15 17

LOT

CX 023

CAL

(Exp.) 2013-03-05
(YYYY - MM - DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	ABX VETPACK										TOLERANCES TOLERANCE
		MICROS	MICROS									
		Scil Vet ABC	Scil Vet ABC Plus ESV 60									
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.4	9.3									± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.45	4.42									± 0.06
HB HGB	g/dl	13.5	13.3									± 0.2
	g/l	135	133									± 2
	mmol/l	8.38	8.26									± 0.12
HT HCT	%	38.7	38.7									± 1.0
	l/l	0.387	0.387									± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	235	236									± 10
VPM MPV	µm ³ ; fl	N/A	8.5									± 0.5

99300778-A
ref: TEMP-0387 Rev.30