

ABX Minocal



LOT CX 013
Rev 1

CAL

(Exp.) 2013-02-05
(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX VETPACK										TOLERANCES TOLERANCE
		MICROS					MICROS					
		Scil Vet ABC	Scil Vet ABC Plus ESV 60									
GB WBC	$10^3/mm^3 ; 10^9/l$	9.2	9.1									± 0.2
GR RBC	$10^6/mm^3 ; 10^{12}/l$	4.54	4.52									± 0.06
HB HGB	g/dl	13.7	13.6									± 0.2
	g/l	137	136									± 2
	mmol/l	8.51	8.45									± 0.12
HT HCT	%	39.4	39.1									± 1.0
	l/l	0.394	0.391									± 0.010
PLAQ. PLTS	$10^3/mm^3 ; 10^9/l$	245	236									± 10
VPM MPV	$\mu m^3 ; fl$	N/A	7.9									± 0.5

Ref: TEMP-0387 Rev.22 FRONT / RECTO 9930044-B