

# ABX Minocal



**LOT** CX 122  
Rev 1

**CAL**

(Exp.) 2013-01-05  
(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX VETPACK										TOLERANCES TOLERANCE
		MICROS					MICROS					
		Scil Vet ABC	Scil Vet ABC Plus ESV 60									
GB WBC	$10^3/mm^3 ; 10^9/l$	<b>9.3</b>	<b>9.3</b>									$\pm 0.2$
GR RBC	$10^9/mm^3 ; 10^{12}/l$	<b>4.58</b>	<b>4.53</b>									$\pm 0.06$
HB HGB	g/dl	<b>13.6</b>	<b>13.6</b>									$\pm 0.2$
	g/l	<b>136</b>	<b>136</b>									$\pm 2$
	mmol/l	<b>8.45</b>	<b>8.45</b>									$\pm 0.12$
HT HCT	%	<b>39.3</b>	<b>39.1</b>									$\pm 1.0$
	l/l	<b>0.393</b>	<b>0.391</b>									$\pm 0.010$
PLAQ. PLTS	$10^3/mm^3 ; 10^9/l$	<b>252</b>	<b>244</b>									$\pm 10$
VPM MPV	$\mu m^3 ; fl$	N/A	<b>7.9</b>									$\pm 0.5$

Ref: TEMP-0387 Rev.22 FRONT / RECTO 9930044-B