

ABX Minocal



LOT CX 102
Rev 1

CAL

(Exp.) 2012-11-05
(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX VETPACK										TOLERANCES TOLERANCE	
		MICROS											
		Scil Vet ABC	Scil Vet ABC Plus ESV 60										
GB WBC	$10^3/mm^3 ; 10^9/l$	9.4	9.6										± 0.2
GR RBC	$10^6/mm^3 ; 10^{12}/l$	4.72	4.72										± 0.06
HB HGB	g/dl	13.6	13.5										± 0.2
	g/l	136	135										± 2
	mmol/l	8.45	8.38										± 0.12
HT HCT	%	39.2	39.4										± 1.0
	l/l	0.392	0.394										± 0.010
PLAQ. PLTS	$10^3/mm^3 ; 10^9/l$	255	247										± 10
VPM MPV	$\mu m^3 ; fl$	N/A	7.7										± 0.5

Ref: TEMP-0387 Rev.22 FRONT / RECTO 9930044-B