

ABX Minocal



LOT CX 024
Rev 2

CAL

 (Exp.) 2014-03-05
(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX VETPACK										TOLERANCES TOLERANCE
		MICROS	MICROS									
		Scil Vet ABC	Scil Vet ABC Plus ESV 60									
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.0	9.1									± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.41	4.41									± 0.06
HB HGB	g/dl	13.5	13.5									± 0.2
	g/l	135	135									± 2
	mmol/l	8.38	8.38									± 0.12
HT HCT	%	38.1	38.4									± 1.0
	l/l	0.381	0.384									± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	250	250									± 10
VPM MPV	µm ³ ; fl	N/A	7.9									± 0.5

Ref: TEMP-0387 Rev.35 9930078-A