

ABX Minocal



LOT CX 072
Rev 1

CAL

(Exp.) 2012-08-05
(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX VETPACK										TOLERANCES TOLERANCE		
		MICROS												
		Scil Vet ABC	Scil Vet ABC Plus ESV 60											
GB WBC	$10^3/mm^3 ; 10^9/l$	9.0	9.1											± 0.2
GR RBC	$10^6/mm^3 ; 10^{12}/l$	4.50	4.52											± 0.06
HB HGB	g/dl	13.5	13.5											± 0.2
	g/l	135	135											± 2
	mmol/l	8.38	8.38											± 0.12
HT HCT	%	38.7	38.9											± 1.0
	l/l	0.387	0.389											± 0.010
PLAQ. PLTS	$10^3/mm^3 ; 10^9/l$	248	250											± 10
VPM MPV	$\mu m^3 ; fl$	N/A	7.8											± 0.5

Ref: TEMP-0387 Rev.22 FRONT / RECTO 9930044-B