

ABX Minocal



LOT CX 052
Rev 1

CAL

(Exp.) 2012-06-05
(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX VETPACK										TOLERANCES TOLERANCE
		MICROS										
		Scil Vet ABC	Scil Vet ABC Plus ESV 60									
GB WBC	$10^3/mm^3 ; 10^9/l$	9.1	9.1									± 0.2
GR RBC	$10^6/mm^3 ; 10^{12}/l$	4.53	4.53									± 0.06
HB HGB	g/dl	13.5	13.3									± 0.2
	g/l	135	133									± 2
	mmol/l	8.38	8.26									± 0.12
HT HCT	%	39.0	38.7									± 1.0
	l/l	0.390	0.387									± 0.010
PLAQ. PLTS	$10^3/mm^3 ; 10^9/l$	250	247									± 10
VPM MPV	$\mu m^3 ; fl$	N/A	7.9									± 0.5

Ref: TEMP-0387 Rev.21 FRONT / RECTO 9930044-B