

# ABX Minocal



**LOT** CX 042  
Rev 1

**CAL**

(Exp.) 2012-05-05  
(YYYY - MM - DD)

		ABX VETPACK										
PARAMETRES PARAMETERS	UNITES UNITS	MICROS	MICROS									TOLERANCES TOLERANCE
		Scil Vet ABC	Scil Vet ABC Plus ESV 60									
GB WBC	10 <sup>3</sup> /mm <sup>3</sup> ; 10 <sup>9</sup> /l	9.0	9.2									± 0.2
GR RBC	10 <sup>6</sup> /mm <sup>3</sup> ; 10 <sup>12</sup> /l	4.49	4.58									± 0.06
HB HGB	g/dl	13.5	13.3									± 0.2
	g/l	135	133									± 2
	mmol/l	8.38	8.28									± 0.12
HT HCT	%	39.5	39.4									± 1.0
	l/l	0.395	0.394									± 0.010
PLAQ. PLTS	10 <sup>9</sup> /mm <sup>3</sup> ; 10 <sup>9</sup> /l	242	240									± 10
VPM MPV	µm <sup>3</sup> ; fl	N/A	7.9									± 0.5

Ref: FORM-0387 Rev.20 FRONT / RECTO 9930044-B