

ABX Minocal



LOT CX 032
Rev 1

CAL

(Exp.) 2012-04-05
(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX VETPACK										TOLERANCES TOLERANCE	
		MICROS											
		Scil Vet ABC	Scil Vet ABC Plus ESV 60										
GB WBC	$10^3/mm^3 ; 10^9/l$	8.8	9.0										± 0.2
GR RBC	$10^6/mm^3 ; 10^{12}/l$	4.51	4.61										± 0.06
HB HGB	g/dl	13.4	13.3										± 0.2
	g/l	134	133										± 2
	mmol/l	8.32	8.28										± 0.12
HT HCT	%	39.0	39.2										± 1.0
	l/l	0.390	0.392										± 0.010
PLAQ. PLTS	$10^3/mm^3 ; 10^9/l$	242	243										± 10
VPM MPV	$\mu m^3 ; fl$	N/A	8.0										± 0.5

Ref: FORM-0387 Rev.20 FRONT / RECTO 9930044-B