

ABX Minocal



LOT CX 121
Rev 1

CAL

(Exp.) 2012-01-05
(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX VETPACK								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC	10 ³ /mm ³ ; 10 ⁹ /l	8.7	8.8							± 0.2
GR RBC	10 ⁹ /mm ³ ; 10 ¹² /l	4.49	4.56							± 0.06
HB HGB	g/dl	13.5	13.4							± 0.2
	g/l	135	134							± 2
	mmol/l	8.38	8.34							± 0.12
HT HCT	%	39.2	39.5							± 1.0
	l/l	0.392	0.395							± 0.010
PLAQ. PLTS	10 ⁹ /mm ³ ; 10 ⁹ /l	260	257							± 10
VPM MPV	µm ³ ; fl	N/A	8.2							± 0.5

Ref: FORM-0387 Rev.20 FRONT / RECTO 9930044-B