

# ABX Minocal



**LOT** CX 034  
Rev 1

**CAL**

 (Exp.) 2014-04-05  
(YYYY - MM - DD)

| PARAMETRES<br>PARAMETERS | UNITES<br>UNITS                                       | ABX VETPACK  |                             |  |  |  |  |  |  |  |  | TOLERANCES<br>TOLERANCE |
|--------------------------|---|--------------|-----------------------------|--|--|--|--|--|--|--|--|-------------------------|
|                          |   | MICROS       | MICROS                      |  |  |  |  |  |  |  |  |                         |
|                          |   | Scil Vet ABC | Scil Vet ABC Plus<br>ESV 60 |  |  |  |  |  |  |  |  |                         |
| GB WBC                   | 10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l  | 9.1          | 9.1                         |  |  |  |  |  |  |  |  | ± 0.2                   |
| GR RBC                   | 10 <sup>6</sup> /mm <sup>3</sup> ;10 <sup>12</sup> /l | 4.66         | 4.62                        |  |  |  |  |  |  |  |  | ± 0.06                  |
| HB HGB                   | g/dl  | 13.6         | 13.5                        |  |  |  |  |  |  |  |  | ± 0.2                   |
|                          | g/l   | 136          | 135                         |  |  |  |  |  |  |  |  | ± 2                     |
|                          | mmol/l  | 8.45         | 8.38                        |  |  |  |  |  |  |  |  | ± 0.12                  |
| HT HCT                   | %   | 39.0         | 38.9                        |  |  |  |  |  |  |  |  | ± 1.0                   |
|                          | l/l   | 0.390        | 0.389                       |  |  |  |  |  |  |  |  | ± 0.010                 |
| PLAQ. PLTS               | 10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l  | 235          | 240                         |  |  |  |  |  |  |  |  | ± 10                    |
| VPM MPV                  | µm <sup>3</sup> ; fl                                  | N/A          | 8.0                         |  |  |  |  |  |  |  |  | ± 0.5                   |

Ref: TEMP-0387 Rev.35 9930078-A