

ABX Minocal

LOT

CX 066

CAL



(Exp.) 2016-07-05
(YYYY - MM - DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio								TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA		PENTRA
			CRP 200			60 8P 60 18P					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l		9.1			9.0					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l		4.52			4.53					± 0.06
HB HGB	g/dl		13.4			13.4					± 0.2
	g/l		134			134					± 2
	mmol/l		8.34			8.31					± 0.12
HT HCT	%		38.4			38.6					± 1.0
	l/l		0.384			0.386					± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l		245			250					± 10
VPM MPV	µm ³ ; fl		8.2			8.2					± 0.5

ABX Minocal

LOT

CX 066

Rev 1

CAL



(Exp.) 2016-07-05
(YYYY-MM-DD)

PARAMETRES PARAMETERS	UNITES UNITS	LyseVet								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC	10 ⁹ /mm ³ ;10 ⁹ /l	8.6	8.8							± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.55	4.58							± 0.06
HB HGB	g/dl	13.6	13.5							± 0.2
	g/l	136	135							± 2
	mmol/l	8.45	8.38							± 0.12
HT HCT	%	38.7	38.7							± 1.0
	l/l	0.387	0.387							± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	237	253							± 10
VPM MPV	µm ³ ; fl	N/A	8.4							± 0.5

Ref: TEMP-0387 Rev.42 FRONT / RECTO 1300019741-A