

# ABX Minocal

**LOT**

**CX 046**

**CAL**



**(Exp.) 2016-05-05**

Rev 1

(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio								TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA		PENTRA
			CRP 200			60 8P 60 18P					
GB WBC	10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l		<b>9.5</b>			<b>9.2</b>					± 0.2
GR RBC	10 <sup>6</sup> /mm <sup>3</sup> ;10 <sup>12</sup> /l		<b>4.70</b>			<b>4.62</b>					± 0.06
HB HGB	g/dl		<b>13.7</b>			<b>13.5</b>					± 0.2
	g/l		<b>137</b>			<b>135</b>					± 2
	mmol/l		<b>8.52</b>			<b>8.37</b>					± 0.12
HT HCT	%		<b>40.0</b>			<b>39.3</b>					± 1.0
	l/l		<b>0.400</b>			<b>0.393</b>					± 0.010
PLAQ. PLTS	10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l		<b>253</b>			<b>240</b>					± 10
VPM MPV	µm <sup>3</sup> ; fl		<b>9.1</b>			<b>8.7</b>					± 0.5

# ABX Minocal

**LOT**

**CX 046**

Rev 1

**CAL**



**(Exp.) 2016-05-05**

(YYYY-MM-DD)

PARAMETRES PARAMETERS	UNITES UNITS	LyseVet								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC	10 <sup>9</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	<b>9.0</b>	<b>9.2</b>							± 0.2
GR RBC	10 <sup>6</sup> /mm <sup>3</sup> ;10 <sup>12</sup> /l	<b>4.70</b>	<b>4.64</b>							± 0.06
HB HGB	g/dl	<b>13.8</b>	<b>13.7</b>							± 0.2
	g/l	<b>138</b>	<b>137</b>							± 2
	mmol/l	<b>8.57</b>	<b>8.51</b>							± 0.12
HT HCT	%	<b>39.7</b>	<b>39.2</b>							± 1.0
	l/l	<b>0.397</b>	<b>0.392</b>							± 0.010
PLAQ. PLTS	10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	<b>240</b>	<b>249</b>							± 10
VPM MPV	µm <sup>3</sup> ; fl	N/A	<b>9.2</b>							± 0.5

Ref: TEMP-0387 Rev.42 FRONT / RECTO 1300019741-A