

ABX Minocal

LOT CX 036
Rev 2

CAL

 (Exp.) 2016-04-05
(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio								TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA		PENTRA
			CRP 200			60 8P 60 18P					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l		9,1			8,9					± 0,2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l		4,46			4,43					± 0,06
HB HGB	g/dl		13,4			13,2					± 0,2
	g/l		134			132					± 2
	mmol/l		8,34			8,19					± 0,12
HT HCT	%		38,6			38,3					± 1,0
	l/l		0,386			0,383					± 0,010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l		250			248					± 10
VPM MPV	µm ³ ; fl		8,6			8,2					± 0,5

ABX Minocal



LOT CX 036
Rev 2

CAL

(Exp.) 2016-04-05
(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	LyseVet								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC	10 ³ /mm ³ ; 10 ⁹ /l	8,7	8,8							± 0,2
GR RBC	10 ⁶ /mm ³ ; 10 ¹² /l	4,44	4,42							± 0,06
HB HGB	g/dl	13,5	13,4							± 0,2
	g/l	135	134							± 2
	mmol/l	8,38	8,32							± 0,12
HT HCT	%	38,6	38,5							± 1,0
	l/l	0,386	0,385							± 0,010
PLAQ. PLTS	10 ³ /mm ³ ; 10 ⁹ /l	240	250							± 10
VPM MPV	µm ³ ; fl	N/A	8,8							± 0,5

Ref: TEMP-0387 Rev.42 FRONT / RECTO 1300019741-A