

ABX Minocal

LOT

CX 026

CAL



(Exp.) 2016-03-05
(YYYY - MM - DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio								TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA		PENTRA
			CRP 200			60 8P 60 18P					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l		9.2			9.1					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l		4.78			4.75					± 0.06
HB HGB	g/dl		13.5			13.4					± 0.2
	g/l		135			134					± 2
	mmol/l		8.40			8.31					± 0.12
HT HCT	%		39.8			39.4					± 1.0
	l/l		0.398			0.394					± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l		258			255					± 10
VPM MPV	µm ³ ; fl		8.3			8.1					± 0.5

ABX Minocal

LOT

CX 026

Rev 1

CAL



(Exp.) **2016-03-05**
(YYYY-MM-DD)

PARAMETRES PARAMETERS	UNITES UNITS	LyseVet								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC	10 ⁹ /mm ³ ;10 ⁹ /l	8.7	8.9							± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.74	4.74							± 0.06
HB HGB	g/dl	13.7	13.6							± 0.2
	g/l	137	136							± 2
	mmol/l	8.51	8.45							± 0.12
HT HCT	%	39.3	39.3							± 1.0
	l/l	0.393	0.393							± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	245	257							± 10
VPM MPV	µm ³ ; fl	N/A	8.5							± 0.5

Ref: TEMP-0387 Rev.40 FRONT / RECTO 1219930095-A