

ABX Minocal

LOT

CX 125

CAL



(Exp.) 2016-01-05
(YYYY - MM - DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio								TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA		PENTRA
			CRP 200			60 8P 60 18P					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l		9.3			9.1					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l		4.46			4.42					± 0.06
HB HGB	g/dl		13.4			13.3					± 0.2
	g/l		134			133					± 2
	mmol/l		8.34			8.25					± 0.12
HT HCT	%		39.0			38.5					± 1.0
	l/l		0.390			0.385					± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l		230			232					± 10
VPM MPV	µm ³ ; fl		8.6			8.0					± 0.5

ABX Minocal

LOT

CX 125

Rev 1

CAL



(Exp.) **2016-01-05**
(YYYY-MM-DD)

PARAMETRES PARAMETERS	UNITES UNITS	LyseVet								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC	10 ⁹ /mm ³ ;10 ⁹ /l	9.1	9.1							± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.50	4.44							± 0.06
HB HGB	g/dl	13.6	13.6							± 0.2
	g/l	136	136							± 2
	mmol/l	8.45	8.45							± 0.12
HT HCT	%	39.1	38.7							± 1.0
	l/l	0.391	0.387							± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	224	230							± 10
VPM MPV	µm ³ ; fl	N/A	8.7							± 0.5

Ref: TEMP-0387 Rev.40 FRONT / RECTO 1219930095-A