

ABX Minocal

LOT

CX 075

CAL



(Exp.) 2015-08-05
(YYYY - MM - DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	Minilysebio								TOLERANCES TOLERANCE	
		MICROS	MICROS	MICROS	MICROS	MICROS	PENTRA	PENTRA	PENTRA		PENTRA
			CRP 200			60 8P 60 18P					
GB WBC	10 ³ /mm ³ ;10 ⁹ /l		9.2			9.0					± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l		4.48			4.49					± 0.06
HB HGB	g/dl		13.6			13.5					± 0.2
	g/l		136			135					± 2
	mmol/l		8.46			8.37					± 0.12
HT HCT	%		39.0			39.0					± 1.0
	l/l		0.390			0.390					± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l		242			240					± 10
VPM MPV	µm ³ ; fl		8.4			7.9					± 0.5

ABX Minocal

LOT

CX 075

Rev 1

CAL



(Exp.) **2015-08-05**
(YYYY-MM-DD)

PARAMETRES PARAMETERS	UNITES UNITS	LyseVet								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC	10 ⁹ /mm ³ ;10 ⁹ /l	8.9	9.0							± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.56	4.52							± 0.06
HB HGB	g/dl	13.8	13.7							± 0.2
	g/l	138	137							± 2
	mmol/l	8.57	8.51							± 0.12
HT HCT	%	39.4	39.1							± 1.0
	l/l	0.394	0.391							± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	232	245							± 10
VPM MPV	µm ³ ; fl	N/A	8.5							± 0.5

Ref: TEMP-0387 Rev.40 FRONT / RECTO 1219930095-A