

ABX Minocal



LOT

CX 055

CAL



(Exp.) 2015-06-05

Rev 1

(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX VETPACK								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.1	9.1							± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.60	4.57							± 0.06
HB HGB	g/dl	13.8	13.8							± 0.2
	g/l	138	138							± 2
	mmol/l	8.57	8.57							± 0.12
HT HCT	%	39.6	39.4							± 1.0
	l/l	0.396	0.394							± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	248	248							± 10
VPM MPV	µm ³ ; fl	N/A	8.6							± 0.5

9930078-A
Ref: TEMP-0387 Rev.38