

ABX Minocal

LOT		CX 045		CAL		(Exp.) 2015-05-05 (YYYY - MM - DD)						
PARAMETRES PARAMETERS	UNITES UNITS	ABX VETPACK									TOLERANCES TOLERANCE	
		MICROS		MICROS								
		Scil Vet ABC	Scil Vet ABC Plus ESV 60									
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.0	9.0								± 0.2	
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.56	4.50								± 0.06	
HB HGB	g/dl	13.6	13.6								± 0.2	
	g/l	136	136								± 2	
	mmol/l	8.45	8.45								± 0.12	
HT HCT	%	39.0	38.7								± 1.0	
	l/l	0.390	0.387								± 0.010	
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	240	242								± 10	
VPM MPV	µm ³ ; fl	N/A	8.7								± 0.5	

 9930078-A
 Ref. TEMP-0387 Rev.38