

ABX Minocal



LOT CX 035
Rev 1

CAL

(Exp.) 2015-04-05
(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX VETPACK								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.1	9.1							± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.58	4.52							± 0.06
HB HGB	g/dl	13.5	13.4							± 0.2
	g/l	135	134							± 2
	mmol/l	8.38	8.32							± 0.12
HT HCT	%	38.7	38.4							± 1.0
	l/l	0.387	0.384							± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	234	240							± 10
VPM MPV	µm ³ ; fl	N/A	8.9							± 0.5

9930078-A
Ref: TEMP-0387 Rev.38