

ABX Minocal



LOT

CX 025

CAL



(Exp.) 2015-03-05

(YYYY - MM - DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	ABX VETPACK								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.4	9.3							± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.50	4.48							± 0.06
HB HGB	g/dl	13.5	13.3							± 0.2
	g/l	135	133							± 2
	mmol/l	8.38	8.26							± 0.12
HT HCT	%	38.4	38.3							± 1.0
	l/l	0.384	0.383							± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	244	245							± 10
VPM MPV	µm ³ ; fl	N/A	8.4							± 0.5

9930078-A
Ref: TEMP-0387 Rev.38