

# ABX Minocal



**LOT**

**CX 015**

**CAL**

 **(Exp.) 2015-02-05**  
(YYYY - MM - DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	ABX VETPACK								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC	10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	<b>9.2</b>	<b>9.3</b>							± 0.2
GR RBC	10 <sup>6</sup> /mm <sup>3</sup> ;10 <sup>12</sup> /l	<b>4.54</b>	<b>4.50</b>							± 0.06
HB HGB	g/dl	<b>13.5</b>	<b>13.3</b>							± 0.2
	g/l	<b>135</b>	<b>133</b>							± 2
	mmol/l	<b>8.38</b>	<b>8.26</b>							± 0.12
HT HCT	%	<b>38.8</b>	<b>38.6</b>							± 1.0
	l/l	<b>0.388</b>	<b>0.386</b>							± 0.010
PLAQ. PLTS	10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	<b>232</b>	<b>242</b>							± 10
VPM MPV	µm <sup>3</sup> ; fl	N/A	<b>8.2</b>							± 0.5

Ref: TEMP-0387 Rev.38 9930078-A