

ABX Minocal



LOT

CX 114

CAL



(Exp.) 2014-12-05

(YYYY - MM - DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	ABX VETPACK								TOLERANCES TOLERANCE
		MICROS								
		Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.2	9.1							± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.50	4.51							± 0.06
HB HGB	g/dl	13.4	13.4							± 0.2
	g/l	134	134							± 2
	mmol/l	8.32	8.32							± 0.12
HT HCT	%	38.5	38.6							± 1.0
	l/l	0.385	0.386							± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	235	242							± 10
VPM MPV	µm ³ ; fl	N/A	9.1							± 0.5

Ref: TEMP-0387 Rev.38 9930078-A