

ABX Minocal



LOT

CX 104

CAL

 (Exp.) **2014-11-05**
(YYYY - MM - DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	ABX VETPACK								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.3	9.3							± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.54	4.58							± 0.06
HB HGB	g/dl	13.6	13.7							± 0.2
	g/l	136	137							± 2
	mmol/l	8.45	8.51							± 0.12
HT HCT	%	38.6	39.0							± 1.0
	l/l	0.386	0.390							± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	235	233							± 10
VPM MPV	µm ³ ; fl	N/A	8.4							± 0.5

Ref: TEMP-0387 Rev.37 9930078-A