

ABX Minocal



LOT

CX 094

CAL



(Exp.) 2014-10-05

(YYYY - MM - DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	ABX VETPACK								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.7	8.8							± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.53	4.54							± 0.06
HB HGB	g/dl	13.6	13.6							± 0.2
	g/l	136	136							± 2
	mmol/l	8.45	8.45							± 0.12
HT HCT	%	38.7	39.0							± 1.0
	l/l	0.387	0.390							± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	235	237							± 10
VPM MPV	µm ³ ; fl	N/A	8.7							± 0.5

Ref: TEMP-0387 Rev.37 9930078-A