

ABX Minocal



LOT CX 084
Rev 1

CAL

 (Exp.) 2014-09-05
(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	ABX VETPACK								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.8	8.9							± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.42	4.45							± 0.06
HB HGB	g/dl	13.5	13.5							± 0.2
	g/l	135	135							± 2
	mmol/l	8.38	8.38							± 0.12
HT HCT	%	38.2	38.7							± 1.0
	l/l	0.382	0.387							± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	248	250							± 10
VPM MPV	µm ³ ; fl	N/A	8.3							± 0.5

Ref: TEMP-0387 Rev.37 9930078-A