

ABX Minocal



LOT

CX 064

CAL



(Exp.) 2014-07-05

(YYYY - MM - DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	ABX VETPACK								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.1	9.2							± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.62	4.60							± 0.06
HB HGB	g/dl	13.8	13.7							± 0.2
	g/l	138	137							± 2
	mmol/l	8.57	8.51							± 0.12
HT HCT	%	39.3	39.3							± 1.0
	l/l	0.393	0.393							± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	244	244							± 10
VPM MPV	µm ³ ; fl	N/A	8.2							± 0.5

Ref: TEMP-0387 Rev.36 9930078-A