

ABX Minocal



LOT

CX 054

CAL



(Exp.) 2014-06-05

(YYYY - MM - DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	ABX VETPACK								TOLERANCES TOLERANCE
		MICROS		MICROS						
		Scil Vet ABC	Scil Vet ABC Plus ESV 60							
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.0	9.1							± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.50	4.48							± 0.06
HB HGB	g/dl	13.7	13.5							± 0.2
	g/l	137	135							± 2
	mmol/l	8.51	8.38							± 0.12
HT HCT	%	38.9	38.8							± 1.0
	l/l	0.389	0.388							± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	246	251							± 10
VPM MPV	µm ³ ; fl	N/A	8.2							± 0.5

Ref: TEMP-0387 Rev.36 9930078-A