

# ABX Minocal



**LOT**

**CX 065**

**CAL**



**(Exp.) 2015-07-05**

(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	WHITEDIFF								TOLERANCES TOLERANCE	
		YUMIZEN									
		H500									
GB WBC	10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	<b>9,0</b>									± 0,2
GR RBC	10 <sup>6</sup> /mm <sup>3</sup> ;10 <sup>12</sup> /l	<b>4,45</b>									± 0,06
HB HGB	g/dl	<b>13,2</b>									± 0,2
	g/l	<b>132</b>									± 2
	mmol/l	<b>8,20</b>									± 0,12
HT HCT	%	<b>37,8</b>									± 1,0
	l/l	<b>0,378</b>									± 0,010
PLAQ. PLTS	10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	<b>259</b>									± 10
VPM MPV	μm <sup>3</sup> ; fl	<b>9,2</b>									± 0,5

Ref: TEMP-0387 Rev.39 BACK / VERSO 1219930095-A