

ABX Minocal



LOT

CX 464

CAL



(Exp.) 2021-12-05

(YYYY - MM - DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	WHITEDIFF								TOLERANCES TOLERANCE
		YUMIZEN	YUMIZEN							
		H550	H500 OT H500 CT							
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.70	8.70							± 0.20
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.51	4.51							± 0.06
HB HGB	g/dl	13.0	13.0							± 0.2
	g/l	130	130							± 2
	mmol/l	8.07	8.07							± 0.12
HT HCT	%	37.4	37.4							± 1.0
	l/l	0.374	0.374							± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	258	258							± 10
VPM MPV	µm ³ ; fl	11.2	11.2							± 0.5

Ref: TEMP-0387 Rev.47 BACK / VERSO 1300060881