

ABX Minocal



LOT

CX 454

CAL



(Exp.) 2021-02-05

(YYYY - MM - DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	WHITEDIFF								TOLERANCES TOLERANCE
		YUMIZEN	YUMIZEN							
		H550	H500 OT H500 CT							
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.75	8.75							± 0.20
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.58	4.58							± 0.06
HB HGB	g/dl	13.3	13.3							± 0.2
	g/l	133	133							± 2
	mmol/l	8.26	8.26							± 0.12
HT HCT	%	38.7	38.7							± 1.0
	l/l	0.387	0.387							± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	255	254							± 10
VPM MPV	µm ³ ; fl	11.3	11.9							± 0.5

Ref: TEMP-0387 Rev.46 BACK / VERSO 1300060881