

# ABX Minocal



**LOT**

**CX 434**

**CAL**



**(Exp.) 2019-06-05**

Rev 1

(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	WHITEDIFF								TOLERANCES TOLERANCE
		YUMIZEN	YUMIZEN							
		H550	H500 OT H500 CT							
GB WBC	10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	<b>8.80</b>	<b>8.90</b>							± 0.20
GR RBC	10 <sup>6</sup> /mm <sup>3</sup> ;10 <sup>12</sup> /l	<b>4.55</b>	<b>4.55</b>							± 0.06
HB HGB	g/dl	<b>13.0</b>	<b>13.0</b>							± 0.2
	g/l	<b>130</b>	<b>130</b>							± 2
	mmol/l	<b>8.07</b>	<b>8.07</b>							± 0.12
HT HCT	%	<b>37.5</b>	<b>37.5</b>							± 1.0
	l/l	<b>0.375</b>	<b>0.375</b>							± 0.010
PLAQ. PLTS	10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	<b>258</b>	<b>262</b>							± 10
VPM MPV	µm <sup>3</sup> ; fl	<b>11.9</b>	<b>10.8</b>							± 0.5

Ref: TEMP-0387 Rev.46 BACK / VERSO 1300060881