

ABX Minocal



LOT

CX 465

CAL



(Exp.) 2022-01-05

Rev 1

(YYYY - MM - DD)

PARAMETRES PARAMETERS	UNITES UNITS	WHITEDIFF								TOLERANCES TOLERANCE
		YUMIZEN	YUMIZEN							
		H550	H500 OT H500 CT							
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	8.80	8.80							± 0.20
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.45	4.45							± 0.06
HB HGB	g/dl	13.1	13.1							± 0.2
	g/l	131	131							± 2
	mmol/l	8.14	8.14							± 0.12
HT HCT	%	37.6	37.6							± 1.0
	l/l	0.376	0.376							± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	267	267							± 10
VPM MPV	µm ³ ; fl	11.4	11.4							± 0.5

Ref: TEMP-0387 Rev.47 BACK / VERSO 1300060881