

ABX Minocal



LOT

CX 414

CAL

(Exp.) **2017-10-05**
(YYYY - MM - DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	WHITEDIFF								TOLERANCES TOLERANCE
		YUMIZEN	YUMIZEN							
		H550	H500 OT H500 CT							
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.0	9.0							± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.47	4.49							± 0.06
HB HGB	g/dl	13.4	13.3							± 0.2
	g/l	134	133							± 2
	mmol/l	8.32	8.26							± 0.12
HT HCT	%	38.5	38.2							± 1.0
	l/l	0.385	0.382							± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	265	270							± 10
VPM MPV	µm ³ ; fl	10.5	10.1							± 0.5

Ref: TEMP-0387 Rev.44 BACK / VERSO 1300032436-A