

ABX Minocal



LOT

CX 404

CAL

(Exp.) **2016-12-05**
(YYYY - MM - DD)

Rev **1**

PARAMETRES PARAMETERS	UNITES UNITS	WHITEDIFF								TOLERANCES TOLERANCE
		YUMIZEN	YUMIZEN							
		H500 OT v1 only	H500 OT H500 CT							
GB WBC	10 ³ /mm ³ ;10 ⁹ /l	9.2	9.3							± 0.2
GR RBC	10 ⁶ /mm ³ ;10 ¹² /l	4.59	4.62							± 0.06
HB HGB	g/dl	13.3	13.4							± 0.2
	g/l	133	134							± 2
	mmol/l	8.26	8.32							± 0.12
HT HCT	%	37.9	38.3							± 1.0
	l/l	0.379	0.383							± 0.010
PLAQ. PLTS	10 ³ /mm ³ ;10 ⁹ /l	242	254							± 10
VPM MPV	µm ³ ; fl	9.6	9.8							± 0.5

Ref: TEMP-0387 Rev.43 BACK / VERSO 1300027265-A