

# ABX Minocal



**LOT**

**CX 463**

**CAL**



**(Exp.) 2021-11-05**

(YYYY - MM - DD)

Rev 1

PARAMETRES PARAMETERS	UNITES UNITS	WHITEDIFF								TOLERANCES TOLERANCE
		YUMIZEN	YUMIZEN							
		H550	H500 OT H500 CT							
GB WBC	10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	<b>8.80</b>	<b>8.80</b>							± 0.20
GR RBC	10 <sup>6</sup> /mm <sup>3</sup> ;10 <sup>12</sup> /l	<b>4.40</b>	<b>4.40</b>							± 0.06
HB HGB	g/dl	<b>12.8</b>	<b>12.8</b>							± 0.2
	g/l	<b>128</b>	<b>128</b>							± 2
	mmol/l	<b>7.95</b>	<b>7.95</b>							± 0.12
HT HCT	%	<b>37.0</b>	<b>37.0</b>							± 1.0
	l/l	<b>0.370</b>	<b>0.370</b>							± 0.010
PLAQ. PLTS	10 <sup>3</sup> /mm <sup>3</sup> ;10 <sup>9</sup> /l	<b>250</b>	<b>250</b>							± 10
VPM MPV	µm <sup>3</sup> ; fl	<b>9.4</b>	<b>9.4</b>							± 0.5

Ref: TEMP-0387 Rev.47 BACK / VERSO 1300060881